

NEUROFILAMENT LIGHT POLYPEPTIDE RABBIT MAB

货号: N261781

产品全名: Neurofilament light polypeptide 兔单克隆抗体

基因符号 NFL; NF-L; NF68; CMT1F; CMT2E; CMTDIG; PPP1R110

UNIPROT ID: P08551

背景: Neurofilaments usually contain three intermediate filament proteins: L, M, and H which are involved in the maintenance of neuronal caliber. Miscellaneous NF-L is the most abundant of the three neurofilament proteins and, like the other nonepithelial intermediate filament proteins, it can form homopolymeric 10-nm filaments.

抗原: Recombinant protein of mouse 68kDa Neurofilament

经过测试的应用: WB, IHC-P, IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-7H3

分子量: Calculated MW: 62 kDa; Observed MW: 68 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Mouse

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience & Mature Neurons

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of 68kDa Neurofilament in mouse brain lysates using Neurofilament light polypeptide antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
