

NEUROKININ 1 RECEPTOR RABBIT PAB

货号: N225949

产品全名: Neurokinin 1 Receptor 兔多抗

基因符号 TACR1; NK1R; TAC1R; Substance-P receptor; SPR; NK-1 receptor; NK-1R; Tachykinin receptor 1

UNIPROT ID: P25103

背景: This is a receptor for the tachykinin neuropeptide substance P. It is probably associated with G proteins that activate a phosphatidylinositol-calcium second messenger system. The rank order of affinity of this receptor to tachykinins is: substance P > substance K > neuromedin-K.

抗原: A synthesized peptide derived from human Neurokinin 1 Receptor

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 46 kDa; Observed MW: 46 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Neurokinin 1 Receptor in K562 lysates using Neurokinin 1 Receptor antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
