

NEUROMEDIN B RECEPTOR RABBIT PAB

货号: N225532

产品全名: Neuromedin B Receptor 兔多抗

基因符号 Neuromedin B receptor antibody; Bombesin Receptor 1 antibody; BBI antibody

UNIPROT ID: P28336

背景: Receptor for neuromedin-B.

抗原: Synthetic Peptide of NMBR

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 43 kDa; Observed MW: 50 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience&Neuromedin B Receptor

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Neuromedin B Receptor in Human Brain, rat Brain lysates using NMBR(Bombesin Receptor)antibody.