

## NFIB RABBIT MAB

货号: N262603

产品全名: NFIB 兔单克隆抗体

基因符号 CTF; HMGIC/NFIB; NF-1/B; NF1-B; NF1B; NFIB2; NFI-B; NFI-RED;  
Nfib; NFIB2; NFIB3

**UNIPROT ID:** O00712

背景: Recognizes and binds the palindromic sequence 5'-TTGGC'NGCCAA-3' present in viral and cellular promoters and in the origin of replication of adenovirus type 2. These proteins are individually capable of activating transcription and replication.

抗原: A synthetic peptide of human NFIB

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R03-0E5

分子量: Calculated MW: 47 kDa; Observed MW: 47,55 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Rat

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of NFIB in C6, HeLa lysates using NFIB antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---