

NICASTRIN RABBIT MAB

货号: N262609

产品全名: Nicastrin 兔单克隆抗体

基因符号 APH2; Ncstn; NCT; Nicastrin

UNIPROT ID: Q92542

背景: Essential subunit of the gamma-secretase complex, an endoprotease complex that catalyzes the intramembrane cleavage of integral membrane proteins such as Notch receptors and APP (beta-amyloid precursor protein). It probably represents a stabilizing cofactor required for the assembly of the gamma-secretase complex.

抗原: A synthetic peptide of human Nicastrin

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-9E6

分子量: Calculated MW: 78 kDa; Observed MW: 110-120 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience&Alzheimer's disease

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Nicastrin in K562, rat Brain, 3T3, HeLa lysates using Nicastrin antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
