

NOXIN RABBIT PAB

货号: N225194

产品全名: NOXIN 兔多抗

基因符号 NOXIN; C11orf82; Nitric oxide-inducible gene 蛋白

UNIPROT ID: Q8IXT1

背景: May be an anti-apoptotic protein involved in DNA repair or cell survival.

抗原: The antiserum was produced against synthesized peptide derived from the Internal region of human NOXIN. AA range:111-160

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 112 kDa; Observed MW: 85 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of NOXIN in 823, 293T, AD293, HeLa, HepG2 lysates using NOXIN antibody.