

NTPCR RABBIT MAB

货号: N262630

产品全名: NTPCR 兔单克隆抗体

基因符号 Clorf57; HCR-NTPase

UNIPROT ID: Q9BSD7

背景: Has nucleotide phosphatase activity towards ATP, GTP, CTP, TTP and UTP. Hydrolyzes nucleoside diphosphates with lower efficiency.

抗原: Recombinant protein of human Nucleoside triphosphate phosphohydrolase

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R09-3II

分子量: Calculated MW: 21 kDa; Observed MW: 21 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Nucleoside triphosphate Phospho-hydrolase in 293T, Jurkat, HeLa lysates using NTPCR antibody.