

## NTT4 RABBIT PAB

货号: N225164

产品全名: NTT4 兔多抗

基因符号 SLC6A17; NTT4; Sodium-dependent neutral amino acid transporter SLC6A17; Sodium-dependent neurotransmitter transporter NTT4; Solute carrier family 6 member 17

**UNIPROT ID:** Q9HIV8

背景: Functions as a sodium-dependent vesicular transporter selective for proline, glycine, leucine and alanine.

抗原: Synthesized peptide derived from NTT4 . at AA range: 240-320

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 81 kDa; Observed MW: 81 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: -

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of NTT4 in rat stomach lysates using NTT4 antibody.