

## NUP107 RABBIT MAB

货号: N263463

产品全名: NUP107 兔单克隆抗体

基因符号 ODG6; NUP84; NPHS11; ODG6; GAMOS7

**UNIPROT ID:** P57740

背景: Plays a role in the nuclear pore complex (NPC) assembly and/or maintenance (PubMed:12552102, PubMed:15229283). Required for the assembly of peripheral proteins into the NPC (PubMed:15229283, PubMed:12552102). May anchor NUP62 to the NPC (PubMed:15229283).

抗原: Recombinant protein of human Nup107

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R06-5F6

分子量: Calculated MW: 106 kDa; Observed MW: 106 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Tags & Cell Markers

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Nup107 in A549, HL-60, U251 lysates using NUP107 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---