

NUP155 RABBIT MAB

货号: N263464

产品全名: NUP155 兔单克隆抗体

基因符号 NI55; ATFB15

UNIPROT ID: O75694

背景: Essential component of nuclear pore complex.

抗原: Recombinant protein of human NUP155

经过测试的应用: WB,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R03-6F7

分子量: Calculated MW: 155 kDa; Observed MW: 155 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Tags & Cell Markers

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of NUP155 in HeLa, A549, HL-60 lysates using NUP155 antibody.



Immunocytochemistry analysis of NUP155 (green) in HeLa using NUP155 antibody, and DAPI (blue).