

NUP50 RABBIT MAB

货号: N263466

产品全名: NUP50 兔单克隆抗体

基因符号 NPAP60; NPAP60L

UNIPROT ID: Q9UKX7

背景: Component of the nuclear pore complex that has a direct role in nuclear protein import.

抗原: A synthetic peptide of human NUP50

经过测试的应用: WB,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-3H1

分子量: Calculated MW: 50 kDa; Observed MW: 50 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: -

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of NUP50 (green) in HeLa using NUP50 antibody, and DAPI (blue)



Western blot analysis of NUP50 in HeLa, C6 lysates using NUP50 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
