

OAT1 RABBIT PAB

货号: N225245

产品全名: OAT1 兔多抗

基因符号 SLC22A6; OAT1; PAHT; Solute carrier family 22 member 6;
Organic anion transporter 1; hOAT1; PAH transporter; hPAHT; Renal organic
anion transporter 1; hROAT1

UNIPROT ID: Q4U2R8

背景: Involved in the renal elimination of endogenous and exogenous
organic anions.

抗原: Synthesized peptide derived from human OAT1 Polyclonal

经过测试的应用: WB, ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 62 kDa; Observed MW: 55 kDa

亚型: IgG

纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.3 containing 50% glycerol,
0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing



Western blot analysis of OAT1 in HEPG2 lysates using OAT1 antibody.