

P2Y12 RABBIT MAB

货号: N263475

产品全名: P2Y12 兔单克隆抗体

基因符号 HORK3; P2Y12; ADPG-R; BDPLT8; SP1999; P2T(AC); P2Y(AC); P2Y(12)R; P2Y(ADP); P2Y(cyc)

UNIPROT ID: Q9H244

背景: Receptor for ADP and ATP coupled to G-proteins that inhibit the adenylyl cyclase second messenger system.

抗原: A synthetic peptide of human P2Y12

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-2E1

分子量: Calculated MW: 39 kDa; Observed MW: 52 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience&P2Y12

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of P2Y12 in A549, HL-60, U251, U87-MG lysates using P2Y12 antibody.