

P73 RABBIT MAB

货号: N262662

产品全名: p73 兔单克隆抗体

基因符号 TP73; P73; Tumor protein p73; p53-like transcription factor; p53-related 蛋白

UNIPROT ID: O15350

背景: p73 IS a protein of the p53 family of proteins. Participates in the apoptotic response to DNA damage. Phosphorylated in a cell cycle-dependent manner and negatively regulated by CDKs. When overproduced, activates transcription from p53-responsive promoters and induces apoptosis. May be a tumor suppressor protein. Seven alternatively-spliced isoforms have been described.

抗原: A synthetic peptide of human p73

经过测试的应用: WB,IP,ChIP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20 ChIP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R07-1J7

分子量: Calculated MW: 70 kDa; Observed MW: 80,70 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of p73 in 293 lysates using p73 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
