

PAOX RABBIT PAB

货号: N225055

产品全名: PAOX 兔多抗

基因符号 PAOX; PAO; Peroxisomal N(1)-acetyl-spermine/spermidine oxidase; Polyamine oxidase

UNIPROT ID: Q6QHF9

背景: Flavoenzyme which catalyzes the oxidation of N1-acetylspermine to spermidine and is thus involved in the polyamine back-conversion.

抗原: Synthesized peptide derived from PAOX . at AA range: 260-340

经过测试的应用: WB,IHC-P,ELISA

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 56 kDa; Observed MW: 70 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction


储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human liver using PAOX antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human liver using PAOX antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of PAOX in rat stomach, K562 lysates using PAOX antibody.