

PAR4 RABBIT MAB

货号: N263480

产品全名: PAR4 兔单克隆抗体

基因符号 PAWR; PAR4; PRKC apoptosis WTI regulator protein; Prostate apoptosis response 4 protein; Par-4

UNIPROT ID: Q96IZ0

背景: Pro-apoptotic protein capable of selectively inducing apoptosis in cancer cells, sensitizing the cells to diverse apoptotic stimuli and causing regression of tumors in animal models. Induces apoptosis in certain cancer cells by activation of the Fas prodeath pathway and coparallel inhibition of NF-kappa-B transcriptional activity.

抗原: A synthetic peptide of human PAR4

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-8C0

分子量: Calculated MW: 37 kDa; Observed MW: 41 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Mouse


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PAR4 in HeLa, A549, HL-60, U251, U87-MG lysates using PAR4 antibody.



Western blot analysis of PAR4 in 3T3 lysates using PAR4 antibody.