

PARATHYROID HORMONE RABBIT PAB

货号: N225837

产品全名: Parathyroid Hormone 兔多抗

基因符号 hPTH; Parathormone; Parathyrin; Parathyroid hormone 1;
Prepro PTH; Preproparathyroid hormone; PTH1; PTH1 receptor; PTHIR; PTHR1

UNIPROT ID: P01270

背景: PTH elevates calcium level by dissolving the salts in bone and preventing their renal excretion. Stimulates [1-14C]-2-deoxy-D-glucose (2DG) transport and glycogen synthesis in osteoblastic cells.

抗原: A synthesized peptide derived from human Parathyroid Hormone

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 13 kDa; Observed MW: 36 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience&Parathyroid Hormone

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Parathyroid Hormone in Human Parathyroid Hormone recombinant protein using Parathyroid Hormone antibody.