

PCBP2 RABBIT MAB

货号: N262678

产品全名: PCBP2 兔单克隆抗体

基因符号 HNRPE2; HNRNPE2; hnRNP-E2

UNIPROT ID: Q15366

背景: Single-stranded nucleic acid binding protein that binds preferentially to oligo dC. Major cellular poly(rC)-binding protein. Binds also poly(rU). Negatively regulates cellular antiviral responses mediated by MAVS signaling (PubMed:19881509). It acts as an adapter between MAVS and the E3 ubiquitin ligase ITCH, therefore triggering MAVS ubiquitination and degradation (PubMed:19881509).

抗原: A synthetic peptide of human PCBP2

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-2F5

分子量: Calculated MW: 39 kDa; Observed MW: 35-45 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Mouse,Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Microbiology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PCBP2/hnRNP E2 in C6, 3T3 lysates using PCBP2 antibody.

