

PCGF6 RABBIT MAB

货号: N262679

产品全名: PCGF6 兔单克隆抗体

基因符号 RING finger protein 134; Mel18 and Bmi1-like RING finger

UNIPROT ID: Q99NA9

背景: Transcriptional repressor. May modulate the levels of histone H3K4Me3 by activating KDM5D histone demethylase. Component of a Polycomb group (PcG) multiprotein PRC1-like complex, a complex class required to maintain the transcriptionally repressive state of many genes, including Hox genes, throughout development.

抗原: Recombinant protein of mouse PCGF6

经过测试的应用: WB, ICC/IF, IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-7J6

分子量: Calculated MW: 40 kDa; Observed MW: 40 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Mouse, Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PCGF6 in PC-12 lysates using PCGF6 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
