

## PCK2 RABBIT PAB

货号: N225907

产品全名: PCK2 兔多抗

基因符号 PCK2; PEPCK2; Phosphoenolpyruvate carboxykinase [GTP]; mitochondrial; PEPCK-M; Phosphoenolpyruvate carboxylase

**UNIPROT ID:** Q16822

**背景:** Catalyzes the conversion of oxaloacetate (OAA) to phosphoenolpyruvate (PEP), the rate-limiting step in the metabolic pathway that produces glucose from lactate and other precursors derived from the citric acid cycle.

**抗原:** A synthesized peptide derived from human PCK2

**经过测试的应用:** WB,IHC-P,ICC/IF,FC

**推荐稀释比:** WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 71 kDa; Observed MW: 71 kDa

**亚型:** IgG


**纯化:** Affinity Chromatography

**种属反应性:** Human, Mouse and Rat

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Cardiovascular

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PCK2 in (1) HepG2 lysates; (2) mouse brain lysates; (3) rat spleen lysates using PCK2 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---