

PDGFR BETA RABBIT MAB

货号: N261763

产品全名: PDGFR beta 兔单克隆抗体

基因符号 PDGFRB; PDGFR; PDGFR1; Platelet-derived growth factor receptor beta; PDGF-R-beta; PDGFR-beta; Beta platelet-derived growth factor receptor; Beta-type platelet-derived growth factor receptor; CD140 antigen-like family member B; Platelet-deri

UNIPROT ID: P09619

背景: This gene encodes a cell surface tyrosine kinase receptor for members of the platelet-derived growth factor family. These growth factors are mitogens for cells of mesenchymal origin. The identity of the growth factor bound to a receptor monomer determines whether the functional receptor is a homodimer or a heterodimer, composed of both platelet-derived growth factor receptor alpha and beta polypeptides.

抗原: A synthetic peptide of human PDGFR beta

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-5B2

分子量: Calculated MW: 124 kDa; Observed MW: 190 kDa

亚型: IgG

纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PDGFR beta in 3T3 lysates using PDGFR beta antibody.



Western blot analysis of PDGFR beta in SH-SY5Y, C6 lysates using PDGFR beta antibody.