

PDIA2 RABBIT PAB

货号: N225774

产品全名: PDIA2 兔多抗

基因符号 P4HB; PDIA2; PDIR

UNIPROT ID: Q13087

背景: Acts as an intracellular estrogen-binding protein. May be involved in modulating cellular levels and biological functions of estrogens in the pancreas. May act as a chaperone that inhibits aggregation of misfolded proteins.

抗原: A synthesized peptide derived from human PDI

经过测试的应用: WB,FC

推荐稀释比: WB: 1/500-1/1000 FC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 58 kDa; Observed MW: 58 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PDI in HeLa lysates using PDIA2 antibody.