

PDLIM1 RABBIT MAB

货号: N262689

产品全名: PDLIM1 兔单克隆抗体

基因符号 PDLIM1; CLIM1; CLP36; PDZ and LIM domain protein 1; C-terminal LIM domain protein 1; Elfin; LIM domain protein CLP-36

UNIPROT ID: O00151

背景: Cytoskeletal protein that may act as an adapter that brings other proteins (like kinases) to the cytoskeleton.

抗原: A synthetic peptide of human PDLIM1

经过测试的应用: WB, ICC/IF, IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-7A8

分子量: Calculated MW: 36 kDa; Observed MW: 36 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of PDLIM1 (green) in HeLa using PDLIM1 antibody, and DAPI (blue)



Western blot analysis of PDLIM1 in K562, C6, 3T3, HeLa lysates using PDLIM1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
