

## PEROXIREDOXIN 4 RABBIT MAB

货号: N262699

产品全名: Peroxiredoxin 4 兔单克隆抗体

基因符号 AOE37-2; PRX-4

**UNIPROT ID:** Q13162

背景: Probably involved in redox regulation of the cell. Regulates the activation of NF-kappa-B in the cytosol by a modulation of I-kappa-B-alpha phosphorylation.

抗原: A synthetic peptide of human Peroxiredoxin 4

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R09-7F9

分子量: Calculated MW: 31 kDa; Observed MW: 27 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Peroxiredoxin 4 in 293T lysates using Peroxiredoxin 4 antibody.