

PEROXIREDOXIN 6 RABBIT MAB

货号: N263490

产品全名: Peroxiredoxin 6 兔单克隆抗体

基因符号 Peroxiredoxin-6; Antioxidant protein 2; Liver 2D page spot 40; NSGPx; PRDX6; AOP2; KIAA

UNIPROT ID: P30041

背景: Involved in redox regulation of the cell. Can reduce H₂O₂ and short chain organic, fatty acid, and phospholipid hydroperoxides. May play a role in the regulation of phospholipid turnover as well as in protection against oxidative injury.

抗原: A synthetic peptide of human Peroxiredoxin 6

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R03-7F5

分子量: Calculated MW: 25 kDa; Observed MW: 25 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Peroxiredoxin 6 in A549, HL-60, U2OS lysates using Peroxiredoxin 6 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
