

## PEX16 RABBIT PAB

货号: N225509

产品全名: PEX16 兔多抗

基因符号 Peroxisomal membrane protein PEX16; Peroxin-16

**UNIPROT ID:** Q9Y5Y5

**背景:** peroxisomal biogenesis factor 16(PEX16) Homo sapiens The protein encoded by this gene is an integral peroxisomal membrane protein. An inactivating nonsense mutation localized to this gene was observed in a patient with Zellweger syndrome of the complementation group CGD/CG9. Expression of this gene product morphologically and biochemically restores the formation of new peroxisomes, suggesting a role in peroxisome organization and biogenesis. Alternative splicing has been observed for this gene and two variants have been described. [provided by RefSeq, Jul 2008],

**抗原:** Synthesized peptide derived from human protein . at AA range: 100-180

**经过测试的应用:** WB,ELISA

**推荐稀释比:** WB: 1/500-1/1000 ELISA: 1/10000

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 39 kDa; Observed MW: 36 kDa

**亚型:** IgG


**纯化:** Affinity Purified

**种属反应性:** Human,Mouse

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** -

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of of PEX16 in PC-12 lysates using PEX16 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---