

## PFKL RABBIT MAB

货号: N262703

产品全名: PFKL 兔单克隆抗体

基因符号 6-phosphofructokinase type B; PFK-B; Phosphohexokinase

**UNIPROT ID:** P17858

背景: Catalyzes the phosphorylation of D-fructose 6-phosphate to fructose 1,6-bisphosphate by ATP, the first committing step of glycolysis.

抗原: A synthetic peptide of human PFKL

经过测试的应用: WB, ICC/IF

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-7A4

分子量: Calculated MW: 85 kDa; Observed MW: 85 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Mouse, Rat

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PFKL in C6, 3T3 lysates using PFKL antibody.