

PHLDA1 RABBIT MAB

货号: N262710

产品全名: PHLDA1 兔单克隆抗体

基因符号 PQ-rich protein; PQR protein; Proline- and histidine-rich protein; T-cell death-associated gene 51 蛋白

UNIPROT ID: Q8WV24

背景: Seems to be involved in regulation of apoptosis.

抗原: A synthetic peptide of human PHLDA1

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-9F4

分子量: Calculated MW: 45 kDa; Observed MW: 45 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PHLDA1 in A375 lysates using PHLDA1 antibody.