

## PHOSPHO-4E BP1 (THR37) RABBIT PAB

货号: N225301

产品全名: Phospho-4E BP1 (Thr37) 兔多抗

基因符号 EIF4EBP1; Eukaryotic translation initiation factor 4E-binding protein 1; 4E-BP1; eIF4E-binding protein 1; Phosphorylated heat- and acid-stable protein regulated by insulin 1; PHAS-I

**UNIPROT ID:** Q13541

**背景:** Translation repressor protein 4E-BP1 (also known as PHAS-1) inhibits cap-dependent translation by binding to the translation initiation factor eIF4E. Hyperphosphorylation of 4E-BP1 disrupts this interaction and results in activation of cap-dependent translation. Both the PI3 kinase/Akt pathway and FRAP/mTOR kinase regulate 4E-BP1 activity.

**抗原:** Synthetic peptide of human EIF4EBP1

**经过测试的应用:** WB,IHC-P

**推荐稀释比:** WB: 1/500-1/1000 IHC: 1/50-1/100

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 13 kDa; Observed MW: 18 kDa

**亚型:** IgG

**纯化:** Affinity Purified


**种属反应性:** Human, Mouse and Rat

**Modification:** Phosphorylated


**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Epigenetics and Nuclear Signaling

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma, using 4EBP1 (Phospho-Thr36) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Western blot analysis of Phospho-4E BP1 (Thr37) in MDAMB435 lysates using Phospho-4E BP1 (Thr37) antibody.