

PHOSPHO-C-JUN (THR93) RABBIT PAB

货号: N226073

产品全名: Phospho-c-Jun (Thr93) 兔多抗

基因符号 Transcription factor AP-1; Activator protein 1; AP1; Proto-oncogene c-Jun; V-jun avian sarcoma virus 17 oncogene homolog; p39

UNIPROT ID: P05412

背景: c-jun is a transcription factor that recognizes and binds to the enhancer heptamer motif 5'-TGA[CG]TCA-3'. Promotes activity of NR5A1 when phosphorylated by HIPK3 leading to increased steroidogenic gene expression upon cAMP signaling pathway stimulation.

抗原: A synthetic phosphopeptide corresponding to residues surrounding Thr93 of human c-Jun

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 36 kDa; Observed MW: 36 kDa

亚型: IgG

纯化: Affinity Purified


种属反应性: Human,Mouse

Modification: Phosphorylated

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Phospho-c-Jun (Thr93) in 3T3, 293, PC-12 lysates using Phospho-c-Jun (Thr93) antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
