

PHOSPHO-C MYB (SER11) RABBIT PAB

货号: N225602

产品全名: Phospho-c Myb (Ser11) 兔多抗

基因符号 MYB; Transcriptional activator Myb; Proto-oncogene c-Myb

UNIPROT ID: P10242

背景: Transcriptional activator; DNA-binding protein that specifically recognize the sequence 5'-YAAC[GT]G-3'. Plays an important role in the control of proliferation and differentiation of hematopoietic progenitor cells.

抗原: A synthesized peptide derived from human Myb

经过测试的应用: WB,IHC-P,ICC/IF,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 72 kDa; Observed MW: 72,80 kDa

亚型: IgG

纯化: Affinity Chromatography

种属反应性: Human

Modification: Phosphorylated

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Phospho-vMyb + c-Myb (S11) in Ramos lysates using Phospho-c Myb (Ser11) antibody.