

PHOSPHO-CSK (SER364) RABBIT PAB

货号: N225328

产品全名: Phospho-CSK (Ser364) 兔多抗

基因符号 CSK; Tyrosine-protein kinase CSK; C-Src kinase; Protein-tyrosine kinase CYL

UNIPROT ID: P41240

背景: Specifically phosphorylates 'Tyr-504' on LCK, which acts as a negative regulatory site. Can also act on the LYN and FYN kinases.

抗原: Synthetic peptide of human CSK

经过测试的应用: WB, ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 51 kDa; Observed MW: 55 kDa

亚型: IgG

纯化: Affinity Purified


种属反应性: Human, Mouse and Rat

Modification: Phosphorylated

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Phospho-CSK (Ser364) in HUVEC lysates using Phospho-CSK (Ser364) antibody.