

PHOSPHO-ERK1/2 (THR202/THR185) RABBIT MAB

货号: N261710

产品全名: Phospho-ERK1/2 (Thr202/Thr185) 兔单克隆抗体

基因符号 MAPK1/MAPK3

UNIPROT ID: P27361/P28482

背景: Serine/threonine kinase which acts as an essential component of the MAP kinase signal transduction pathway. MAPK1/ERK2 and MAPK3/ERK1 are the 2 MAPKs which play an important role in the MAPK/ERK cascade. They participate also in a signaling cascade initiated by activated KIT and KITLG/SCF. Depending on the cellular context, the MAPK/ERK cascade mediates diverse biological functions such as cell growth, adhesion, survival and differentiation through the regulation of transcription, translation, cytoskeletal rearrangements.

抗原: A synthetic phosphopeptide corresponding to residues surrounding Thr185 of human ERK2

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R05-2H9

分子量: Calculated MW: 44,42 kDa; Observed MW: 44,42 kDa

亚型: IgG

纯化: Affinity Purified


种属反应性: Human,Hamster,Rat

Modification: Phosphorylated


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Phospho-p44/42 MAPK (ERK1/2) (Thr202) in C6, CHO-K1 lysates using Phospho-p44/42 MAPK (ERK1/2) (Thr202) antibody.



Western blot analysis of Phospho-ERK1 (Thr202)/ERK2 (Thr185) in HeLa, HeLa+TPA lysates using Phospho-ERK1/2 (Thr202/Thr185) antibody.