

PHOSPHO-JAK1 (TYR1022) RABBIT PAB

货号: N225383

产品全名: Phospho-JAK1 (Tyr1022) 兔多抗

基因符号 JAK1; JAK1A; JAK1B; Tyrosine-protein kinase JAK1; Janus kinase 1; JAK-1

UNIPROT ID: P23458

背景: Tyrosine kinase of the non-receptor type, involved in the IFN-alpha/beta/gamma signal pathway. Kinase partner for the interleukin (IL)-2 receptor.

抗原: The antiserum was produced against synthesized peptide derived from human JAK1 around the phosphorylation site of Tyr1022. AA range:988-1037

经过测试的应用: ICC/IF, WB, IHC-F, IHC-P, ELISA

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 133 kDa; Observed MW: 133 kDa

亚型: IgG

纯化: Affinity Purified


种属反应性: Human, Mouse and Rat

Modification: Phosphorylated

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Phospho-JAK1 (Tyr1022) in various lysates using Phospho-JAK1 (Tyr1022) antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
