

## PHOSPHO-PAK4/5/6 (SER474/SER560/SER602) RABBIT PAB

货号: N225909

产品全名: Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) 兔多抗

基因符号 PAK4; PAK5; PAK6

**UNIPROT ID:** O96013/Q9P286/Q9NQ5

背景: Serine/threonine protein kinase that plays a role in a variety of different signaling pathways including cytoskeleton regulation, cell migration, growth, proliferation or cell survival.

抗原: A synthesized peptide derived from human Phospho-PAK4 + PAK5 + PAK6 (S474 + S560 + S602)

经过测试的应用: WB,FC

推荐稀释比: WB: 1/500-1/1000 FC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 64,80,75 kDa; Observed MW: 64,80,75 kDa

亚型: IgG

纯化: Affinity Chromatography


种属反应性: Human, Mouse and Rat

**Modification:** Phosphorylated

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Phospho-PAK4 + PAK5 + PAK6 (S474 + S560 + S602) in U87-MG lysates using Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---