

PHOSPHO-PKA RII ALPHA (SER99) RABBIT PAB

货号: N225559

产品全名: Phospho-PKA RII alpha (Ser99) 兔多抗

基因符号 PRKAR2A; PKR2; PRKAR2; cAMP-dependent protein kinase type II-alpha regulatory subunit

UNIPROT ID: P13861

背景: Regulatory subunit of the cAMP-dependent protein kinases involved in cAMP signaling in cells. Type II regulatory chains mediate membrane association by binding to anchoring proteins, including the MAP2 kinase.

抗原: A synthesized peptide derived from human Phospho-PKA R2 (S99)

经过测试的应用: WB,IHC-P,ICC/IF,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 46 kDa; Observed MW: 51 kDa

亚型: IgG

纯化: Affinity Chromatography


种属反应性: Human, Mouse and Rat,Pig

Modification: Phosphorylated

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Phospho-PKA R2 (Ser99) in K562 lysates using Phospho-PKA RII alpha (Ser99) antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
