

## PHOSPHO-VIMENTIN (SER56) RABBIT MAB

货号: N263070

产品全名: Phospho-Vimentin (Ser56) 兔单克隆抗体

基因符号 VIM; Vimentin

**UNIPROT ID:** P08670

**背景:** Vimentin an intermediate filament protein. Intermediate filament proteins are expressed in a tissue-specific manner. Desmin is the subunit specific for muscle and vimentin the subunit specific for mesenchymal tissue.

**抗原:** A synthetic phosphopeptide corresponding to residues surrounding Ser56 of human Vimentin

**经过测试的应用:** WB, ICC/IF, IP

**推荐稀释比:** WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

**种属反应性:** Rabbit

**克隆性:** Rabbit Monoclonal

**克隆编号:** R01-0E2

**分子量:** Calculated MW: 54 kDa; Observed MW: 54 kDa

**亚型:** IgG

**纯化:** Affinity Purified


**种属反应性:** Human

**Modification:** Phosphorylated


**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Neuroscience

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Phospho-Vimentin (Ser56) in U251MG lysates using Phospho-Vimentin (Ser56) antibody.



Immunocytochemistry analysis of Phospho-Vimentin (Ser56) (green) in HeLa using Phospho-Vimentin (Ser56) antibody, and DAPI (blue).