

PHOSPHO-VIMENTIN (SER72) RABBIT MAB

货号: N263071

产品全名: Phospho-Vimentin (Ser72) 兔单克隆抗体

基因符号 VIM; Vimentin

UNIPROT ID: P08670

背景: Vimentin an intermediate filament protein. Intermediate filament proteins are expressed in a tissue-specific manner. Desmin is the subunit specific for muscle and vimentin the subunit specific for mesenchymal tissue.

抗原: A synthetic phosphopeptide corresponding to residues surrounding Ser72 of human Vimentin

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-1C8

分子量: Calculated MW: 54 kDa; Observed MW: 54 kDa

亚型: IgG

纯化: Affinity Purified


种属反应性: Rat

Modification: Phosphorylated

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Vimentin (Phospho-Ser72) in rat Brain lysates using Phospho-Vimentin (Ser72) antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
