

## PHOSPHO-VIMENTIN (SER72) RABBIT PAB

货号: N225562

产品全名: Phospho-Vimentin (Ser72) 兔多抗

基因符号 VIM; Vimentin

**UNIPROT ID:** P08670

**背景:** Vimentin an intermediate filament protein. Intermediate filament proteins are expressed in a tissue-specific manner. Desmin is the subunit specific for muscle and vimentin the subunit specific for mesenchymal tissue.

**抗原:** A synthesized peptide derived from human Phospho-Vimentin (S72)

**经过测试的应用:** WB,IP

**推荐稀释比:** WB: 1/500-1/1000 IP: 1/20

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 54 kDa; Observed MW: 54 kDa

**亚型:** IgG

**纯化:** Affinity Chromatography


**种属反应性:** Human, Mouse and Rat

**Modification:** Phosphorylated

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Neuroscience

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Phospho-Vimentin (Ser72) in HeLa lysates treated with Calyculin A using Phospho-Vimentin (Ser72) antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---