

PI3 KINASE C2 GAMMA RABBIT PAB

货号: N225203

产品全名: PI3 Kinase C2 gamma 兔多抗

基因符号 PIK3C2G; Phosphatidylinositol 4-phosphate 3-kinase C2 domain-containing subunit gamma; PI3K-C2-gamma; PtdIns-3-kinase C2 subunit gamma; Phosphoinositide 3-kinase-C2-gamma

UNIPROT ID: O75747

背景: Generates phosphatidylinositol 3-phosphate (PtdIns3P) and phosphatidylinositol 3,4-bisphosphate (PtdIns(3,4)P₂) that act as second messengers (By similarity).

抗原: Synthesized peptide derived from the N-terminal region of human PI 3-Kinase C2 γ .

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 166 kDa; Observed MW: 160 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PI3 Kinase C2 gamma in A549 lysates using PI3 Kinase C2 gamma antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
