

## PI3 KINASE P110 BETA RABBIT PAB

货号: N225668

产品全名: PI3 Kinase p110 beta 兔多抗

基因符号 PIK3CB; DKFZp779K1237; MGC133043; PI3K; PI3KCB; PI3Kbeta; PIK3C1; p110-BETA

**UNIPROT ID:** P42338

**背景:** Phosphoinositide 3-kinase (PI3K) catalyzes the production of phosphatidylinositol-3,4,5-triphosphate by phosphorylating phosphatidylinositol (PI), phosphatidylinositol-4-phosphate (PIP) and phosphatidylinositol-4,5-bisphosphate (PIP2). Growth factors and hormones trigger this phosphorylation event, which in turn coordinates cell growth, cell cycle entry, cell migration, and cell survival.

**抗原:** A synthesized peptide derived from human PI3 Kinase p110 beta

**经过测试的应用:** WB,IP,FC

**推荐稀释比:** WB: 1/500-1/1000 IP: 1/20 FC: 1/50-1/100

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 123 kDa; Observed MW: 110 kDa

**亚型:** IgG


**纯化:** Affinity Chromatography

**种属反应性:** Human

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Cell Biology

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PI3 Kinase p110 beta in Jurkat lysates using PI3 Kinase p110 beta antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---