

PIAS1/2 RABBIT MAB

货号: N263618

产品全名: PIAS1/2 兔单克隆抗体

基因符号 PIAS1; PIAS2; PIAS3; Protein inhibitor of activated STAT protein 1

UNIPROT ID: O75925/O75928/Q9Y6X2

背景: O75928,Q9Y6X2.This proteins are involved in the pathway protein sumoylation, which is part of Protein modification.

抗原: A synthesized peptide derived from human PIAS1

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R06-7C1

分子量: Calculated MW: 72 kDa; Observed MW: 76 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PIAS1+PIAS2 in HeLa, A549, HL-60, U2OS lysates using PIAS1+PIAS2 antibody.



Western blot analysis of PIAS1/PIAS2 in mouse testis, rat testis lysates using PIAS1/2 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
