

PIM2 RABBIT MAB

货号: N263494

产品全名: PIM2 兔单克隆抗体

基因符号 PIM2; Serine/threonine-protein kinase pim-2; Pim-2h

UNIPROT ID: Q9PIW9

背景: Promotes cell survival in response to a variety of proliferative signals via positive regulation of the I-kappaB kinase/NF-kappaB cascade; this process requires phosphorylation of MAP3K8/COT. Prevents apoptosis induced by growth factor withdrawal via inhibition of caspase-3 activation, and via phosphorylation of pro-apoptotic proteins.

抗原: A synthetic peptide of human PIM2

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-5D8

分子量: Calculated MW: 34 kDa; Observed MW: 40,38,34 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PIM2 in A549, HL-60, U251 lysates using PIM2 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
