

## PKA C GAMMA RABBIT MAB

货号: N262758

产品全名: PKA C gamma 兔单克隆抗体

基因符号 PRKACG; cAMP-dependent protein kinase catalytic subunit gamma; PKA C-gamma

**UNIPROT ID:** P22612

背景: Phosphorylates a large number of substrates in the cytoplasm and the nucleus.

抗原: Recombinant protein of human PRKACG

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-8II

分子量: Calculated MW: 40 kDa; Observed MW: 40 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PRKACG in MCF-7, C6, HeLa lysates using PRKACG antibody.



Western blot analysis of PRKACG in mouse testis, mouse lung lysates using PKA C gamma antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---