

PLACENTAL LACTOGEN RABBIT PAB

货号: N225211

产品全名: Placental Lactogen 兔多抗

基因符号 chorionic somatomammotropin hormone 1; PL; CSA; CS-1; CSMT; hCS-1; hCS-A

UNIPROT ID: P0DML2

背景: Produced only during pregnancy and is involved in stimulating lactation, fetal growth and metabolism.

抗原: Synthetic peptide from human protein at AA range: 161-210

经过测试的应用: WB,IHC-P,ELISA

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 25 kDa; Observed MW: 25 kDa

亚型: IgG

纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Placental Lactogen in mouse brain lysates using Placental Lactogen antibody.