

PLASMINOGEN RECEPTOR RABBIT PAB

货号: N225210

产品全名: Plasminogen Receptor 兔多抗

基因符号 Plasminogen receptor (KT; Plg-R(KT))

UNIPROT ID: Q9HBL7

背景: Receptor for plasminogen.

抗原: The antiserum was produced against synthesized peptide derived from the Internal region of human PLGRKT. AA range:1-50

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 17 kDa; Observed MW: 17 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: -

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Plasminogen Receptor in HEPG2 lysates using Plasminogen Receptor antibody.