

## PNKP RABBIT MAB

货号: N263498

产品全名: PNKP 兔单克隆抗体

基因符号 PNK; AOA4; MCSZ; CMT2B2; EIEE10

**UNIPROT ID:** Q96T60

背景: Plays a key role in the repair of DNA damage, functioning as part of both the non-homologous end-joining (NHEJ) and base excision repair (BER) pathways.

抗原: A synthetic peptide of human PNK/PNKP

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R03-0B8

分子量: Calculated MW: 57 kDa; Observed MW: 57 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: -

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using PNK/PNKP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of PNK/PNKP in HeLa, A549, HL-60, U251 lysates using PNKP antibody.