

POLIOVIRUS RECEPTOR RABBIT MAB

货号: N263501

产品全名: Poliovirus Receptor 兔单克隆抗体

基因符号 PVR; PVS; Poliovirus receptor; Nectin-like protein 5; NECL-5; CD155

UNIPROT ID: P15151

背景: Mediates NK cell adhesion and triggers NK cell effector functions. Binds two different NK cell receptors: CD96 and CD226. These interactions accumulates at the cell-cell contact site, leading to the formation of a mature immunological synapse between NK cell and target cell.

抗原: Recombinant protein of human Poliovirus Receptor/PVR

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-2C7

分子量: Calculated MW: 45 kDa; Observed MW: 60-80 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Poliovirus Receptor/PVR antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Poliovirus Receptor/PVR in Hela, A549, U2OS lysates using Poliovirus Receptor antibody.