

## PPCS RABBIT MAB

货号: N263504

产品全名: PPCS 兔单克隆抗体

基因符号 CMD2C

**UNIPROT ID:** Q9HAB8

背景: Catalyzes the second step in the biosynthesis of coenzyme A from vitamin B5, where cysteine is conjugated to 4'-phosphopantothenate to form 4-phosphopantothenoylcysteine.

抗原: A synthetic peptide of human PPCS

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-6I3

分子量: Calculated MW: 34 kDa; Observed MW: 34 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PPCS in HeLa, Jurkat lysates using PPCS antibody.