

PPP2R5D RABBIT MAB

货号: N262749

产品全名: PPP2R5D 兔单克隆抗体

基因符号 B56D; MRD35; B56delta

UNIPROT ID: Q14738

背景: The B regulatory subunit might modulate substrate selectivity and catalytic activity, and also might direct the localization of the catalytic enzyme to a particular subcellular compartment.

抗原: Recombinant protein of human PPP2R5D

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-2A4

分子量: Calculated MW: 70 kDa; Observed MW: 70,66 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PPP2R5D in K562, rat Brain, C6, 3T3, HeLa lysates using PPP2R5D antibody.